






**SC-5 SimCube® NIBP Simulator**



**SC-4 SimCube® NIBP Simulator**



**SC-3 SimCube® NIBP Simulator**

			
	<b>SC-5 SimCube® NIBP Simulator</b>	<b>SC-4 SimCube® NIBP Simulator</b>	<b>SC-3 SimCube® NIBP Simulator</b>
<b>NIBP Modes</b>	Standard, High, Low, Neonatal NIBP Simulation	Standard, High, Neonatal NIBP Simulation	Standard, High, Neonatal NIBP Simulation
<b>Digital Manometer</b>	- 400 to + 400 mmHg	0 to +400 mmHg	0 to +400 mmHg
<b>Leak Test Mode</b>	Yes	No	No
<b>Over Pressure (Peak Detect) Mode</b>	Yes	Yes	Yes
<b>ECG 12 Lead</b>	Yes	Yes	No
<b>Arrhythmia</b>	Yes	No	No
<b>HR Sequence (Alarm Test)</b>	Yes	Yes	No
<b>IBP</b>	Yes	No	No
<b>Optional Features</b>	Battery Boost Option, Serial Communications for expanded NIBP simulations and automation	Battery Boost Option	Battery Boost Option