

FlowTrax FT-2 Infusion Pump Analyzer Features:

Unbeatable Value

- The fastest, most accurate infusion pump analyzer on the market
- NEW Quick Prime — faster priming with no stopping to re-prime between tests!
- First Multifunction IV Pump Analyzer: eliminate expense and calibration of four pieces of test equipment with one device
- Field Serviceable: the only analyzer with removable sample tube and input luer for easy maintenance.
- All features are NIST traceable.
- Ultra-durable: passes drop testing 50X from 3 ft!
- Operates on two 'AA' batteries
- 4-year warranty

UNSURPASSED Rate and Volume Accuracy

- Results within 1% accuracy for rates from 1 ml/hr to 999 ml/hr in less than 3 minutes
- NEW Electronic Occlusion feature
- Test any type of infusion device from PCA, syringe, volumetric to feeding pumps
- Unrivalled speed — requires as little as 0.5 ml volume for low rate tests to achieve 1% accuracy — up to 40X faster than other products
- Never miss a test result with Trend feature
- Capture real-time results to any File, Form or Record on your PC with DataSnap
- Market-best measurement resolution of 0.0008 milliliters

Full-Featured Pressure and Vacuum Meter

- 0.5% Accuracy of Reading
- Utilize as standalone pressure meter. Test any type of medical device for pressure or vacuum (-12.5 to +75 PSI).
- NEW Electronic Occlusion Feature (0-40 PSI) with graphical chart and total elapsed time of test
- Measure real time and peak pressure
- 9 units of measure: mmHg, PSI (gauge), inH₂O, cmH₂O, kPa, atm, inHg, BAR (absolute and gauge)

Temperature Meter

- 0.1% ± 0.2 °C Accuracy
- Compatible with any YSI 400 temperature probe
- Measure fluid, gas or air temperatures (-28.9 to +115.5 °C)
- Test incubators, fluid warmers, freezers and more

Digital Timer / Stop watch

- Integrated digital timer and stop watch
- Automatically starts and stops with fluid measurements.