

**SC-5 SimCube® NIBP Simulator Features:**

- Adult (120/80 mmHg), Neonatal (70/40 mmHg), Hypertensive (190/120 mmHg), and Hypotensive (80/40 mmHg) NIBP Simulations
- Digital Manometer
- 12-lead ECG Simulation
- Arrhythmia
- Respiration
- Invasive Blood Pressure (IBP) Dynamic values (120/80, 70/40, 190/120)
- 7 static IBP values (0, 25, 50, 100, 150, 200, 250)
- Heart Rate/Respiration sequence
- Pacer
- Over Pressure (Peak Detect) Test mode measures the peak pressure produced by the patient monitor. The readout latches the highest pressure value sensed, in 0.1mmHg steps.
- Leak Rate Test Mode

**Package Contents:**

- (1) SC-5 SimCube® NIBP Simulator
- (1) Power Supply (6VDC)
- (1) Certificate of Calibration
- (1) NIBP Adapters

Optional Accessories: NIBP Adapters, Carrying Case, Cuff Jacket, ECG Extender Module, Battery Boost, IBP adapters